

# Power Scrub Vacuum Line Cleaner Directions for Use

Solmetex  
**POWERSCRUB™**  
VACUUM LINE CLEANER



**Extend your NXT Hg5 warranty from 2 to 5 years—just use PowerScrub, the only vacuum line cleaner we recommend for the Hg5.**

Solmetex recommends daily cleaning of the vacuum lines in each of your operatories.

1. Invert PowerScrub bottle to fill measuring reservoir located at top of container.
2. Tilt container in the opposite direction and drain excess until you have the desired volume of concentrate. Pour into the PowerScrub Sidekick Dispenser as directed below:
  - **For 1 Operatory:**  
Measure 1 fluid ounce of concentrate and mix with 32 oz (1000 ml) of water.
  - **For 2 Operatories:**  
Measure 2 fluid ounces of concentrate and mix with 64 oz (2000 ml) of water.
3. Attach both the chair-side low-volume and high-volume evacuators to the nozzles on the Sidekick Dispenser, REMOVE OR LOOSEN CAP TO VENT, then aspirate 32 oz (1000 ml) of mixed cleaning solution per operatory. There is no need to flush lines after cleaning. Any cleaning solutions left in the dispenser can be poured down the cuspidor.

**First Time Users:** Use half the amount of concentrate or 0.5 fl oz mixed with 32 oz (1000 ml) of water per operatory for the first 2 cleanings to gradually remove accumulated debris from vacuum lines. On subsequent cleanings, follow “full strength” mixture directions above.

**PowerScrub Vacuum Line Cleaner is:**

- ✓ Non-foaming
- ✓ Compatible with all amalgam separator and vacuums.
- ✓ EPA-compliant neutral pH formulation.
- ✓ The only vacuum line cleaner recommended for the NXT Hg5 Amalgam Separator



**Scan here for PowerScrub™ downloads and important information**



**WARNING** ⚠ Keep out of reach of children. Harmful if taken internally. Do not get in eyes, on skin or clothing. Wear eye protection, gloves and protective clothing.

**FIRST AID:** In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. If solution is swallowed, call a physician or poison control immediately. Refer to SDS for complete safety information.

Safe for the environment and septic systems. Non-flammable. Keep from freezing.