



INSTRUCTIONS FOR USE

- Place entire kit into freezer overnight before collecting the samples.
- Scan the order form QR code on the front of the package and complete your office information.

If unable to scan QR code, visit my-solutions-center.solmetex.com/dentisafetesting and CAREFULLY ENTER THE TEST ID NUMBER above the QR code on the front of the package, then continue completing your office information. This alerts the lab to let them know your water sample(s) are being sent.

- Cover workspace with patient napkin and apply clean gloves.
- Flush DUWL for approximately 20 seconds before filling the vial.
- Remove collection vial from kit and place the white cap face down on the napkin.
- Fill collection vial, return white cap and tighten. It is recommended to use a single vial per line tested.
- Add each vial ID code to your digital order form, being careful to input each code accurately.
- **IMPORTANT:** Place frozen icepack into kit, remove adhesive strip and seal the mailer.
- Print return shipping label (available in My Solutions Center, in the Water Testing Tab), attach to the mailer and schedule drop off to UPS.

NOTE: Please DO NOT ship samples on Fridays, weekends or holidays. Only ship test samples on Mondays-Thursdays in order to ensure that samples are viable when received by the lab.

- Results will be posted on the My Solutions Center portal at Solmetex.com. You will receive an email informing you when your lab results are ready to view.



POWERED BY

DENTISAFE
PREMIUM WATER TEST

CUSTOMER SERVICE
1-800-216-5505