


# We have amalgam separation down to a science.

Solmetex has quickly become the industry leader in Amalgam Separation. What makes Solmetex different is our “Total Solution Approach” — integrating the simple design of the NXT Hg5 with waste handling and recycling into our complete product line, providing a truly Green set of solutions for the dental practice.

## Solutions

The NXT Hg5 provides the complete solution:

- New compact design, internal manifold, ideal for hard to fit areas
- Functional with all wet and dry vacuums
- Environmentally friendly packaging
- Simple mail-back recycling program
- Meets all “EPA Dental Rule” requirements
- ISO 11143 certified by NSF International



### Easy

- Clear system, easy to inspect
- Two-minute container change out
- No tools necessary

### Flexible

- Versatile, wall or floor mounted
- Compatible with wet or dry vacuums
- 3 different sizes to fit your practice

### Support

- Compliance center
- Online certificates
- Auto-ship program

### Recycle

- Ship container via UPS
- EPA certified mercury facility
- Recycle Kit includes container and recycle packaging

# Ordering Information

To order your Solmetex NXT Hg5 System, please contact your local dental product dealer. If you don't have a dealer, visit [www.solmetex.com](http://www.solmetex.com) Click on "dealer locator" to find the dental product dealer closest to you or contact Solmetex for information at **1-800-216-5505**.



**NXT Hg5 mini  
Amalgam Separator**

- Item # **NXT-HG5-MINI**  
Specifications:
- 25”h x 11”w x 8”d
  - 3/4” inlet and outlet piping
  - Fernco® couplers
  - Floor / Wall mount bracket
  - Services 1 - 4 chairs
  - Container Reorder # **NXT-HG5-002CR**



**NXT Hg5  
Amalgam Separator**

- Item # **NXT-HG5-001**  
Specifications:
- 30”h x 11”w x 8”d
  - 1-1/2” inlet and outlet piping
  - Fernco® couplers
  - Floor / Wall mount bracket
  - Services 1 - 10 chairs
  - Container Reorder # **NXT-HG5-002CR**



**Hg5 High Volume  
Amalgam Separator**

- Item # **HG5HV-001**  
Specifications:
- 28”h x 18”w x 22.5”d
  - 3” inlet and outlet piping
  - Wall mount bracket
  - Services 11 - 20 chairs
  - Container Reorder # **HG5-002CR**

All Collection Container with Recycle Kits contain:

- Collection Container
- Proper packaging for shipping
- Recycling documentation
- Recycling to EPA approved facility



**NXT Hg5  
Collection Container**

with Recycle Kit  
Item # **NXT-HG5-002CR**



**Hg5  
Collection Container**

with Recycle Kit [legacy]  
Item # **HG5-002CR**



**RAMVAC Utility  
Collection Container**

with Recycle Kit  
Item # **HG5-CR-RVC**

**Solmetex – clearly the leader in amalgam separation.**

The Hg5 amalgam separator is the industry standard and has won more awards than all competitors combined.

