

# Amalgam Waste Recycling Solmetex — The Complete Solution

Market Leader

Clear by Design


ISO 11143 Certified

Easy Collection  
Container Change-Out

**Wet Waste**  
A certified installed NXT Hg5 amalgam separator will capture all the waste and solids sent through your dental water lines before it reaches your vacuum pump. This collection could potentially add 5 years\* to the life of your vacuum.



**NXT Hg5  
Amalgam Separator**



**PWS  
Amalgam  
Waste  
Bucket**

Affordable, One-Stop  
Mail-Back Service

Ensures Compliance  
with New EPA  
Dental Amalgam Rule

No Contracts Required

**Dry Waste**  
A Practice Waste Solutions amalgam bucket will store contact and non-contact amalgam, spent amalgam capsules, chair-side traps, central vacuum filters and extracted teeth containing amalgam. When full, simply print out a shipping label on Solmetex.com and mail it for proper recycling.

**= Compliance**

Certificates of Compliance available  
online 24/7 at Solmetex.com

 **Solmetex**  
Clearly the leader in amalgam separation