

Dealer Catalog Description Page

Hg5 Series of Products (June 2018)

Name: NXT Hg5 Amalgam Separator

- Suitable for 1 to 10 Ops
- New compact design for easier & flexible installation
- Clear by Design
- ISO 11143:2008 certified at 99% separation rate
- New internal manifold design
- New eco-friendly packaging, 100% recyclable
- Item # **NXT-HG5-001**

Name: NXT Hg5 mini Amalgam Separator

- Suitable for 1 to 4 Ops
- New compact design for easier & flexible installation
- Clear by Design
- ISO 11143:2008 certified at 99% separation rate
- New internal manifold design
- New eco-friendly packaging, 100% recyclable
- Item # **NXT-HG5-MINI**

Name: NXT Hg5 High Volume Amalgam Separator

- Suitable for 11 to 20 Ops
- Clear by Design
- ISO 11143 certified
- Fits most large practice/clinic applications
- Item # **NXT-HG5-HV**

Name: NXT Hg5 Collection Container with Recycle Kit

- Replacement NXT Hg5 Collection Container
- Compatible only with the NXT Hg5 Amalgam Separator, the NXT Hg5 mini Amalgam Separator and the NXT Hg5 High Volume Amalgam Separator
- New eco-friendly packaging, 100% recyclable
- New easier shipping procedure
- Same simple mail-back recycling program
- Item # **NXT-HG5-002CR**

Name: Hg5 Collection Container with Recycle Kit [Legacy]

- Replacement Hg5 Collection Container [Legacy]
- Compatible only with the Hg5 Amalgam Separator, the Hg5 mini Amalgam Separator and the Hg5 High Volume Amalgam Separator [All Legacy]
- Recycle packaging
- Simple mail-back recycling program
- Item # **HG5-002CR**

Name: RAMVAC Utility Hg5 Collection Container with Recycle Kit

- Replacement RAMVAC Utility Hg5 Collection Container
- Compatible with the RAMVAC Utility Hg5 Amalgam Separator only
- Recycle packaging
- Simple mail-back recycling program
- Item # **HG5-CR-RVC**

Name: NXT Hg5 Compliance Kit

- Contains one NXT Hg5 Collection Container with Recycle Kit
- Contains one PWS 1.25 gallon Amalgam Bucket
- Collection Container compatible with the NXT Hg5 Amalgam Separator only
- Recycle packaging
- Simple mail-back recycling program
- Item # **NXT-HG5-CK**

Name: Hg5 Compliance Kit

- Contains one Hg5 Collection Container with Recycle Kit
- Contains one PWS 1.25 gallon Amalgam Bucket
- Collection Container compatible with the Hg5 Amalgam Separator only
- Recycle packaging
- Simple mail-back recycling program
- Item # **HG5-CK**

Practice Waste Solutions by Solmetex

June 2018

Name: Solmetex Amalgam Bucket

The Solmetex Amalgam Bucket meets all BMP requirements for the new EPA amalgam regulation. Includes amalgam bucket, recycling, shipping, and recycling certificates.

Items to be recycled in the amalgam bucket

- Contact amalgam/ non-contact amalgam
- Removed teeth with amalgam fillings
- Expended amalgam capsules
- Chairside traps
- Vacuum pump filters

Item # PWS-AB-1	1.25 gallon amalgam bucket
Item # PWS-AB-2	2.5 gallon amalgam bucket
Item # PWS-AB-3.5	3.5 gallon amalgam bucket
Item # PWS-AB-5	5 gallon amalgam bucket

Name: Solmetex Lead Bucket

The Solmetex Lead Bucket meets all BMP requirements for disposing dental lead waste. Includes lead bucket, recycling, shipping, and recycling certificates.

Items to be recycled in lead bucket

- Bitewing x-ray foils
- Sterilization tape
- Lead apron (PWS-A)

Item # PWS- LB-1	1.25 gallon lead bucket
Item # PWS- A	Lead Apron bucket (Maximum of 2)

Name: Solmetex Biohazard/Sharps Solutions

Solmetex Biohazard/Sharps Solutions meets all OSHA requirements for Biohazard and Sharps disposal. Includes container, recycling, shipping, and recycling certificates.

Items placed in biohazard sharps containers

- Sharps Containers*
- Sponges
- Gloves
- Towels
- Gauze
- Bibs

*20 Gallon system can accommodate the following:				
1 Quart up to 28	1 Gallon up to 8	5 Quart up to 6	2 Gallon up to 5	3 Gallon up to 3

Item # **PWS-BH-5** 5 gallon bio-hazard & sharps container disposal
Item # **PWS-BH- 20** *20 gallon bio-hazard & sharps container disposal

Name: Solmetex Photo-Chemical Solutions

Solmetex Photo-Chemical Solutions meets all EPA requirements for Fixer disposal. Includes container, recycling, shipping, and recycling certificates.

Features:

- Easy Installation
- Meets all EPA Requirements for Fixer Disposal

Item # **PWS-PCF-A** Photo Chemical Filter, Automatic W/ Wall Mount
Item # **PWS- PCF-REPL** Photo Chemical Filter, Automatic Replacement Filter
Item # **PWS-PCF-M** Photo Chemical Filter, Manual
Item # **PWS-PCF-CA** California Photo Chemical Filter, Automatic
Item # **PWS-PCF-CM** California Photo Chemical Filter, Manual