

The only  
vacuum line cleaner  
that contains  
**MEREDUCE\***  
Scientifically proven to  
inhibit the release of  
harmful mercury

# Solmetex™ **POWERSCRUB™** VACUUM LINE CLEANER

## Clean Better Flow Better Practice Better...

- Use PowerScrub and extend the warranty of your current NXT Hg5 Amalgam Separator from 2 to 5 years
- SmartScrub Technology — Continuous cleaning action
- Microbial formulation optimizes collection container capacity
- Cleans vacuum line debris on contact
- Releases microbial cleaning enzymes at any water temperature
- Meets all EPA “Dental Rule” requirements — Neutral pH > 6 and < 8
- Compatible with both dry and wet ring vacuums
- Compatible with all amalgam separators



- \* Testing of mercury laden water treated for seven days with MEREDUCE resulted in removal of 98% of mercuric ions.  
— Long Island Analytical Laboratories (Sept 2021)
- \*\* Testing of mercury laden water control samples treated for five days with MEREDUCE resulted in no detection of mercuric ions.  
— Univ of Texas at San Antonio, School of Dentistry (Oct 2021)



# POWERSCRUB™

## VACUUM LINE CLEANER

Solmetex  
recommends daily  
cleaning of your  
vacuum lines

## Ordering Information

Ordering is easy! Simply call Solmetex or place your order through your favorite dental dealer. Questions? Call Solmetex Customer Care 800-216-5505 or chat with us online at [solmetex.com](http://solmetex.com)



**PowerScrub  
Vacuum Line  
Cleaner**

Item # **PCS-VLCB**  
One 64 oz. bottle of  
PowerScrub  
Vacuum Line Cleaner



**PowerScrub  
Vacuum Line Cleaner  
Intro Kit**

Item # **PCS-VLIK**  
• One (64 oz.) bottle of PowerScrub  
Vacuum Line Cleaner  
• One PowerScrub Sidekick Dispenser



**PowerScrub  
Vacuum Line Cleaner  
Sidekick**

Item # **PCS-VLSK**  
The PowerScrub  
Vacuum Line Cleaner  
Sidekick Dispenser



**PowerScrub  
Vacuum Line Cleaner  
Case Kit**

Item # **PCS-VLCK**  
Four (64 oz.) bottles  
of PowerScrub  
Vacuum Line Cleaner

**PowerScrub™ is the only vacuum line  
cleaner recommended for NXT Hg5  
Amalgam Separators**

\* If you are not fully satisfied with the PowerScrub Vacuum Line Cleaner, please contact Solmetex Customer Service at 800-216-5505