



Hg5[®] Series of Amalgam Separators

We've got amalgam separation down to a science. Removing mercury from the dental practice (and the environment) is a good thing - and we're way ahead of the game. When it comes to Amalgam Separation, no one knows more than Solmetex. 21 years ago, our experience in water science led us to create a system that easily and affordably removes amalgam waste from dental practices. Our signature product the Hg5 is the industry standard and has won more awards than all competitors combined. As a 'total solution provider,' Solmetex takes care of waste handling and recycling so you don't have to. Good for you and good for the environment. Solmetex – clearly the leader in amalgam separation

SolmeteX provides the complete solution:

- ISO 11143 certified at greater than 99%
- Professionally engineered, robust, simple design.
- Easy to operate and maintenance free.
- Flexible design for hard-to-fit applications.
- Functional for wet or dry vacuum systems (install before pump on wet vacuum system, install before tank on dry vacuum system).



Product Awards

Dental Advisor recipient of:

Editors' Choice with a 5+ rating
2014 Preferred Product
2015 Preferred Product

Townie Choice Award best in category for Amalgam Separators twelve consecutive years.

dental
townie
choice
awards.2014



1.800.216.5505

www.solmetex.com

Hg5-001

Hg5 Amalgam Separator

Specification Standard:

- 29"h x 10"w x 8"d.
- 1-1/2" inlet and outlet piping.
- Fernco® couplers.
- Floor / Wall mount bracket.
- Services 1 - 10 chairs.

Hg5HV-001

Hg5 High Volume Amalgam Separator

Specification High Volume System:

- 18.5"h x 21"w x 22.5"d.
- 3" inlet and outlet piping.
- Wall mount bracket.
- Services 11 - 20 chairs.

Hg5-Mini

Hg5 Mini Amalam Separator

Specification Mini System:

- 21.5"h x 11"w x 8"d.
- 3/4" inlet and outlet piping.
- Fernco® couplers.
- Floor / Wall mount bracket.
- Services 1 - 4 chairs.

Hg5-002CR

Collection Container with Recycle Kit

Complete container & recycle kit:

- Collection container
- Proper packaging for shipping
- Recycling
- Recycling documentation
- One size fits all systems



To order your SolmeteX Hg5 System, please contact your local dental product dealer. If you don't have a dealer, contact SolmeteX for information at 1-800-216-5505 or visit <http://www.solmetex.com> Click on "dealer locator" to find the dental product dealer closest to you.



Flexible -

Free standing or wall mount
Wet or dry vacuum systems
Multiple sizes to fit your practice

Easy -

Maintenance free
Collection container changed with no tools
Two-minute container change out
Clear collection container lets you see when full

Recycle -

Ship container via UPS to our contracted EPA certified mercury facility
Recycle Kit includes container and recycle packaging plus pre-paid shipping
Certificate provided on-line after recycling

Support -

Expert, responsive, service oriented