

The only
vacuum line cleaner
that contains
MEREDUCE*
Scientifically proven to
inhibit the release of
harmful mercury

Solmetex™ **POWERSCRUB™** VACUUM LINE CLEANER

Clean Better Flow Better Practice Better...

- Use PowerScrub and extend the warranty of your current NXT Hg5 Amalgam Separator from 2 to 5 years
- SmartScrub Technology — Continuous cleaning action, erodes vacuum line debris on contact
- Microbial formulation optimizes collection container capacity
- Releases microbial cleaning enzymes at any water temperature
- Meets all EPA “Dental Rule” requirements — Neutral pH > 6 and < 8
- Compatible with both dry and wet ring vacuums
- Compatible with all amalgam separators



- * Testing of mercury laden water treated for seven days with MEREDUCE resulted in transformation of 98% of mercuric ions.
— Long Island Analytical Laboratories (Sept 2021)
- ** Testing of mercury laden water control samples treated for five days with MEREDUCE resulted in no detection of mercuric ions.
— Univ of Texas at San Antonio, School of Dentistry (Oct 2021)



POWERSCRUB™

VACUUM LINE CLEANER

Solmetex
recommends daily
cleaning of your
vacuum lines

Ordering Information

Ordering is easy! Simply call Solmetex or place your order through your favorite dental dealer. Questions? Call Solmetex Customer Care 800-216-5505 or chat with us online at solmetex.com



**PowerScrub
Vacuum Line
Cleaner**

Item # **PCS-VLCB**
One 64 oz. bottle of
PowerScrub
Vacuum Line Cleaner



**PowerScrub
Vacuum Line Cleaner
Intro Kit**

Item # **PCS-VLIK**
• One (64 oz.) bottle of PowerScrub
Vacuum Line Cleaner
• One PowerScrub Sidekick Dispenser



**PowerScrub
Vacuum Line Cleaner
Sidekick**

Item # **PCS-VLSK**
The PowerScrub
Vacuum Line Cleaner
Sidekick Dispenser



**PowerScrub
Vacuum Line Cleaner
Case Kit**

Item # **PCS-VLCK**
Four (64 oz.) bottles
of PowerScrub
Vacuum Line Cleaner

**PowerScrub™ is the only vacuum line
cleaner recommended for NXT Hg5
Amalgam Separators**