



# STERISIL®

*Safe, Simple, Sterisil*

○ **SHOCK**  
○ **TEST**  
✓ **MAINTAIN**

## Antimicrobial Bottle



Constructed with silver impregnated plastic, the Sterisil® Antimicrobial Bottle helps suppress biofilm development on the internal bottle surface.

### FEATURES AND BENEFITS

- Suppresses the growth of bacteria that may cause odor or deterioration of the bottle itself
- Expanded bottle mouth enables easy removal and replacement of straws
- Molded volumetric indicators for easy water measurement
- Available in 1L and 2L sizes





### THE BACTERIA PROBLEM

Biofilm development inside your dental bottle can contribute to an increase in bacterial content at the dental handpiece.

Visit us online at [Sterisil.com](http://Sterisil.com)




The bottle's silver impregnated construction helps suppress the biofilm development that causes unpleasant odors and/or deterioration of the bottle itself.\*



Use Citrisil™ Shock tablets to clean the bottle's inner surface every time you shock your dental unit waterlines.

Citrisil™ Shock tablets are Sterisil approved and recommended for use with the Sterisil® Antimicrobial Bottle.



Sterisil, Inc.  
835 S. Hwy 105 Ste. D  
Palmer Lake, CO USA 80133  
+1 719 622 7200

*\*The Sterisil® Antimicrobial Bottle is not intended for use as a stand alone water treatment solution.*



For more information contact your **Sterisil® sales representative or local dealer.**  
719 622 7200 or online at [Sterisil.com](http://Sterisil.com)

