



**STERISIL<sup>®</sup>** *Safe, Simple, Sterisil*



# Sterisil<sup>®</sup> Ac System

## EASY, IN-OFFICE AUTOCLAVE WATER PRODUCTION FOR DENTISTRY



### KEY FEATURES

- Eliminates the need for purchasing and storing distilled water
- Output water is suitable for any instrument washer or autoclave
- Recommended for offices with <150ppm TDS in their source water
- Live water quality readings
- Audible and visual alarms
- Autoclave filler wand included
- Self-closing valve heads make cartridge changes mess-free
- Mounts under a cabinet or on the wall

The Sterisil<sup>®</sup> Ac System is designed to deliver purified autoclave water, made right from your municipal supply. No more going to the supermarket or storing giant jugs to keep your autoclaves supplied with water. A simple plug and play solution to one of dentistry's persistent problems.

Visit us online at **Sterisil.com**





STERISIL®  
Safe, Simple, Sterisil

## STOP STORING WATER

Most practices need distilled water and lots of it. But purchasing and storing distilled water is a hassle. Why not make your own? With the Ac, you have an endless supply of distilled quality water without having to procure, store, and handle those large jugs.




## DESIGNED FOR DENTISTRY

With dental specific features like the filler wand, water quality monitoring, and audible and visual alarms, the Ac System delivers where consumer grade units fall short. The Ac System is recommended for offices with  $\leq 150$ ppm TDS content in their tap water. The unit is capable of dispensing purified water containing  $<10$ ppm TDS, making it suitable for use with any autoclave.

STERISIL PART #	DESCRIPTION
<b>AC</b>	Complete Ac System produces autoclave water for areas with water containing $<150$ ppm TDS
<b>AC-C</b>	Single replacement cartridge—option for limited space
<b>AC-CXL</b>	Single replacement XL cartridge
<b>AC-CK</b>	Cartridge replacement kit for Ac System—used for limited space, includes 3 cartridges
<b>AC-CKXL</b>	XL cartridge replacement kit for Ac System, includes 3 cartridges



Proudly made  
in the USA,

 Sterisil, Inc.  
835 S. Hwy 105 Ste. D  
Palmer Lake, CO USA 80133  
+1 719 622 7200



For more information contact your **Sterisil® sales representative or local dealer.**  
719 622 7200 or online at [Sterisil.com](http://Sterisil.com)



GP6000.C