



STERISIL[®] *Safe, Simple, Sterisil*



Sterisil[®] Ac+ System

CONVENIENT, ON DEMAND AUTOCLAVE WATER,
MADE RIGHT IN YOUR PRACTICE



*Tuttnauer Tiva 2H and
T-Edge Autoclaves*

DESIGNED FOR DENTISTRY

- Eliminates the need for purchasing and storing distilled water
- Produces up to 50 gallons of distilled quality water per day
- Output water contains <10ppm TDS and is suitable for any instrument washer or autoclave

ULTRA SAFE, ULTRAVIOLET EQUIPPED

- Class B UV treatment
- Reduces the elevated microbial load that inevitably occurs over time when storing water in tanks
- Helps ensure processed water is safe for dental unit waterlines

BUILT FOR CONVENIENCE

- Can be mounted under a sink or on the wall in a sterilization center
- User serviceable without a dental technician
- Audible and visual maintenance notifications
- Water can be dispensed via a faucet or the Sterisil[®] Wand depending on the practice needs

Visit us online at **Sterisil.com**





STOP BUYING JUGS OF WATER

Most practices need distilled water and lots of it. But purchasing and storing distilled water is a hassle. Why not make your own?

With the Sterisil Ac+, you have an endless supply of distilled quality water without having to procure, store, and handle those large jugs.

The Ac+ connects directly to your municipal water supply and produces up to 50 gallons per day of distilled quality water for your instrument washers, autoclaves, and dental bottles.



ULTRAVIOLET PROTECTION

Due to the slow rate of water production, most dental water systems rely on storage tanks to dispense water on-demand. This stored water inevitably becomes contaminated with bacteria over time. Unless treated, this purified water would be unsuitable for use in dental units.

The Sterisil® Ac+ features a powerful Class B UV treatment to help ensure dispensed water meets the <500 CFU/ml drinking water standard.

Start filling your bottles with confidence, not bacteria-contaminated water.


CONSISTENTLY PURE

Municipal water quality is highly variable depending on your location. This chemical variability can undermine the effectiveness of microbial treatment products like straws and tablets.

The Ac+ removes municipal contaminants and helps facilitate more consistent waterline treatment outcomes.



Proudly made
in the USA.

 Sterisil, Inc.
835 S. Hwy 105 Ste. D
Palmer Lake, CO USA 80133
+1 719 622 7200



For more information contact your **Sterisil® sales representative or local dealer.**

719 622 7200 or online at Sterisil.com

GP7000.C

