



## **SOLMETEX** **COMPLIANCE** **PROGRAM™**

Set your  
schedule  
and we do  
the rest

**Alleviate the guesswork around EPA compliance  
and protect your office from unnecessary fines.  
Enroll in the Solmetex Compliance Program.**

Our new program includes auto-delivery of all the products required by the EPA's "Dental Rule" regulation on amalgam waste. You set the custom schedule that works best for your practice and Solmetex handles the rest.  
*No Contracts. No Fees.*



**COMPLIANCE  
CENTER**



**AMALGAM  
BUCKET**



**COLLECTION  
CONTAINER**



**VACUUM LINE  
CLEANER**

**(877) 207-1551 | [www.solmetex.com](http://www.solmetex.com)**

# Solmetex™ Proudly Introduces the



## Enroll today, be compliant tomorrow

In response to the EPA Dental Rule, Solmetex is excited to introduce its proprietary **Solmetex Compliance Program**. A program designed as a turnkey solution to alleviate guesswork associated with the regulation to ensure your practice is always in compliance. Enrollment enables real-time access to all your practice's documentation and recycling activity online 24/7/365, access to documents that a local inspector may request at any time.

The Solmetex Compliance Program covers the **four main components of the EPA Dental Rule** that offices are required to have to maintain compliance:

### 1. Hg5 Series of Collection Containers with Recycle Kit



**NXT Hg5  
Collection Container**

with Recycle Kit

Item # **NXT-HG5-002CR**

- The only collection container approved for use with the NXT Hg5 Amalgam Separator
- Clear by design for easy visual inspection
- Simple change-out process when solids reach the full line or at 12 months, whichever occurs first.
- Packaging and shipping approved by the US DOT.
- Certificates of Recycling housed and available 24/7 in the Solmetex Compliance Center



**Hg5 [Legacy]  
Collection Container**

with Recycle Kit

Item # **HG5-002CR**

### 2. PowerScrub™ Vacuum Line Cleaner



#### PowerScrub™ Vacuum Line Cleaner Case Kit

Item #

**PCS-VLCK**

- Per EPA regulation requirements, formulation consists of a neutral pH level between 6 and 8.
- Case kit includes four 64oz bottles.
- Recommended for use with the NXT Hg5 Amalgam Separator
- Maximizes capacity of the NXT Hg5 collection container.
- Compatible with all competing amalgam separators
- Compatible with all dry/wet vacuums
- Works at all water temperatures

### 3. Practice Waste Solutions Amalgam Bucket



#### PWS Amalgam Waste Buckets

1.25 Gallon Bucket	Item # <b>PWS-AB-1</b>
2.5 Gallon Bucket	Item # <b>PWS-AB-2</b>
3.5 Gallon Bucket	Item # <b>PWS-AB-3.5</b>
5 Gallon Bucket	Item # <b>PWS-AB-5</b>

- Designed to address the "dry waste" requirements of the EPA Dental Rule
- Per EPA regulation requirements, clearly labeled for amalgam waste.
- Available in four (4) sizes
- Cost of PWS amalgam bucket includes cost of shipping and recycling.
- Certificates of recycling available 24/7 in the Solmetex Compliance Center

### 4. Solmetex Compliance Center



- Available 24/7/365
- Recycling History
- Certificates of Recycling
- Certificates of Installation
- Certificates of Warranty

Participation in the Solmetex Compliance Program is accompanied by a **non-contractual**, automatic fulfillment program for all the products and technologies necessary to ensure compliance with the EPA's regulation on the proper handling of amalgam waste within the dental practice.