

Solmetex Photo-Chemical Solutions

Silver waste from x-ray film processing equipment is filtered, rendering the waste non hazardous to meet local, state, and federal requirements for discharge into public sewer. Spent filters are shipped for recycling via DOT/EPA approved storage/shipping containers. Containers have proper labeling and shipping documents to transport waste for recycling.

- Manual filters
- Automatic filters
- California specific (manual and automatic)
- Easy installation



Suggested Retail

PWS-PCF-A	Photo-Chemical filter, automatic	\$266.00
PWS-PCF-A-REPL	Photo-Chemical filter, replacement	\$140.00
PWS-PCF-M	Photo-Chemical, manual	\$266.00
PWS-PCF-CA	Photo-Chemical, California automatic	\$475.00
PWS-PCF-CM	Photo-Chemical filter, California manual	\$443.00

Biohazard / Sharps Solutions

Collect sharps containers, towels, gloves, articles of clothing, bandages and gauze. Containers have proper labeling and shipping documents to transport waste for recycling.

Items placed in biohazard sharps containers:

- Sharps containers
- Gauze
- Bibs
- Gloves
- Sponges
- Towels



Suggested Retail


PWS-BH-5	5 Gallon container disposal	\$168.00
PWS-BH-20	20 Gallon container disposal*	\$390.00

*20 Gallon system can accommodate the following:

1 Quart up to 28	1 Gallon up to 8	5 Quart up to 6	2 Gallon up to 5	3 Gallon up to 3
---------------------	---------------------	--------------------	---------------------	---------------------



Waste management must be a part of every practice's daily routine. Having a step by step plan in place will ensure a safe and secured environment not only for the dental office but for patients as well.




Our Promise

Solmetex wants to be your "Total Solution Provider" for your dental practice. With our proven experience and end-to-end service, we can help practices run safely with a clean conscience. Our customer support team is ready to assist with the knowledge, the tools and the help. We provide amalgam, biohazard/sharps, lead and photo-chemical products to assist with the steps to dispose and recycle dental waste.

The NXT Hg5 amalgam separator program is an all-in-one solution, assisting the dental office to achieve a regulatory compliant practice with amalgam waste recycling and management.

The Solmetex NXT Hg5 system has received the Townie Choice Award FIFTEEN consecutive years in the Amalgam Separator category, making the Hg5 Series the most award winning amalgam separator in the North America.



 Solmetex™

1-800-216-5505
50 Bearfoot Road • Northborough, MA • 01532
www.solmetex.com



All services and products of Solmetex are subject to change. Trademarks are the property of their respective owners.

Tri-fold 04.2018

 Solmetex™

Clearly The Leader in
Amalgam Separation &
Practice Waste Solutions



Practice Waste Solutions

The Practice Waste Solutions line of products provides your practice the opportunity to meet all your waste disposal and recycling needs from one source. Solmetex offers a variety of disposal kits for amalgam, lead, photo-chemical and biohazard/ sharps in a variety of sizes to suit your practice. These products, offered through your dental dealer, include a simple mail-back recycling program along with certificates of recycling and compliance accessible 24/7 at www.solmetex.com.



Solmetex Amalgam Bucket

Meets all requirements of BMP's and "EPA Dental Rule". Waste amalgam is collected and placed into DOT/EPA approved storage/shipping containers. Containers have proper labeling and shipping documents to transport waste for recycling.

- Chairside traps
- Contact amalgam
- Non-contact amalgam
- Removed teeth with amalgam fillings
- Expended amalgam capsules
- Vacuum pump filters

Suggested Retail

PWS-AB-1	1.25 Gallon amalgam bucket	\$115.00
PWS-AB-2	2.5 Gallon amalgam bucket	\$168.00
PWS-AB-3.5	3.5 Gallon amalgam bucket	\$181.00
PWS-AB-5	5 Gallon amalgam bucket	\$199.00



Solmetex Lead Bucket

Lead waste and lead foil are collected and placed into DOT/EPA approved storage/shipping containers for proper recycling. Containers have proper labeling and shipping documents to transport waste for recycling.

- Bitewing x-ray
- Sterilization tape
- Lead apron bucket (max of 2)

Suggested Retail

PWS-LB-1	1.25 Gallon lead bucket	\$93.00
PWS-LB-A	2.5 Gallon lead bucket (2 aprons max)	\$101.00

NXT Hg5 Series of Amalgam Separators



NXT Hg5 mini (1-4 Chairs) | NXT Hg5 Standard (1-10 Chairs) | NXT Hg5 High Volume (1-20 Chairs)

The Hg5 series meets all requirements of BMP's and "EPA Dental Rule", plus adds value to the equipment room by extending the life of your vacuum system.

The NXT Hg5 provides the complete solution:

- New compact design, internal manifold, ideal for hard to fit areas
- Easy to operate & virtually maintenance-free
- Environmentally friendly packaging, 100% recyclable
- Affordable & cost-efficient
- Functional for wet or dry vacuum systems
- Hg5 series is clear by design
- Supported by responsive experts
- ISO 11143:2008 certified by NSF International
- Simple mail-back recycling program

Suggested Retail

NXT-Hg5-001	NXT Hg5 amalgam separator (container reorder # NXT-Hg5-002CR)	\$860.00
NXT-Hg5-MINI	NXT Hg5 mini amalgam separator (container reorder # NXT-Hg5-002CR)	\$860.00
NXT-Hg5-HV	NXT Hg5 High volume amalgam separator (container reorder # NXT-Hg5-002CR)	\$2,860.00

Backplates included with each system

- NO contracts with hidden fees
- NO pumps
- NO electricity
- NO timers
- NO daily maintenance
- NO additional charges for shipping to recycling facility
- NO decanting

Hg5 Replacement Container

Solmetex Hg5 Recycle Kit mail in program includes a new container, proper box for shipping spent container, shipping to the recycler, recycling and documentation.



Suggested Retail

NXT-Hg5-002CR	NXT Hg5 Collection Container with recycle kit (To be used with the NXT Hg5 Amalgam Separator Only)	\$338.00
Hg5-002CR	Hg5 Collection Container with recycle kit (To be used with the Hg5 Legacy and Hg5 High Volume Amalgam Separators only)	\$338.00
Hg5-CR-RVC	RAMVAC Utility Collection Container with recycle kit (To be used with the RAMVAC Utility Hg5 Amalgam Separator only)	\$338.00

Online certificates

Now you can get all recycling certificates for all products online 24/7 at www.Solmetex.com.



Auto-ship program

Never worry about chair downtime again! We work with you to create a custom schedule to deliver products as you need them.

Contact Solmetex at **800-216-5505** to receive more details and sign up for the Auto-ship program.