



The NXT Hg5 Amalgam Separator

Dentistry's Best Selling, Longest Lasting
Amalgam Separator

Trust the Leader

For over 25 years, Solmetex has led the amalgam separation industry with its award winning NXT Hg5 Amalgam Separator which easily and affordably removes amalgam waste from dental practices.

As a "Total Solution Provider", Solmetex not only takes care of all amalgam waste handling but also meets all regulatory requirements. Plus the NXT Hg5 series of Amalgam Separators add value to the equipment room by extending the life of the vacuum system and is virtually maintenance free.

Solutions

The NXT Hg5 provides the complete solution:

- Available in 3 sizes to suit your practice (mini 1-4 ops; standard 1-10 ops; high volume 11-20 ops)
- Functional with all wet and dry vacuums
- Environmentally friendly packaging
- ISO 11143:2008 certified by NSF International



Flexibility

- Versatile, wall or floor mount
- Compatible with wet or dry vacuums
- May be installed left or right orientation

Ease

- Clear system, easy to inspect
- Two-minute container change out
- No tools necessary

Support

- Full Service Customer Support Department

Ordering Information

There are several ways to order your Solmetex NXT Hg5 System, either call Solmetex **877-207-1551**, visit **www.solmetex.com** to place an online order or contact your local dental product dealer. If you do not have a dealer, please visit solmetex.com and click on "dealer locator" to find a dental dealer near you.

NXT Hg5 Series of Amalgam Separators



**NXT Hg5 mini
Amalgam Separator**

Item # **NXT-HG5-MINI-I**
 Specifications:
 • 25" h x 11" w x 8" d
 • 3/4" inlet & outlet piping
 • Fernco® couplers
 • Floor / Wall mount bracket
 • Services 1 - 4 chairs
 • Container Reorder No.
NXT-HG5-002CRI



**NXT Hg5
Amalgam Separator**

Item # **NXT-HG5-001I**
 Specifications:
 • 30" h x 11" w x 8" d
 • 1 1/2" inlet & outlet piping
 • Fernco® couplers
 • Floor / Wall mount bracket
 • Services 1 - 10 chairs
 • Container Reorder No.
NXT-HG5-002CRI



**NXT Hg5 High Volume
Amalgam Separator**

Item # **NXT-HV-001I**
 Specifications:
 • 29 3/4" h x 17 5/8" w x 22 1/2" d
 • 3" inlet & outlet piping
 • Wall mount bracket
 • Services 11 - 20 chairs
 • Container Reorder No.
NXT-HG5-002CRI

Hg5 Series of Collection Containers



**NXT Hg5
Collection Container**

Item # **NXT-HG5-002I**



**Hg5 [Legacy]
Collection Container**

Item # **HG5-002I**



**RAMVAC Utility
Collection Container**

Item # **HG5CAN-002RVC**

Solmetex -- Clearly the leader in amalgam separation. The NXT Hg5 amalgam separator is the industry standard and has won more awards than all competitors combined.

