

Ordering Information

There are several ways to order your Solmetex NXT Hg5 System, either call Solmetex **877-207-1551**, visit **www.solmetex.com** to place an online order or contact your local dental product dealer. If you do not have a dealer, please visit solmetex.com and click on "dealer locator" to find a dental dealer near you.

NXT Hg5 Series of Amalgam Separators



**NXT Hg5 mini
Amalgam Separator**

Item # **NXT-HG5-MINI**

Specifications:

- 25”h x 11”w x 8”d
- ¾” inlet & outlet piping
- Fernco® couplers
- Floor / Wall mount bracket
- Services 1 - 4 chairs
- Container Reorder No. **NXT-HG5-002CR**



**NXT Hg5
Amalgam Separator**

Item # **NXT-HG5-001**

Specifications:

- 30”h x 11”w x 8”d
- 1½” inlet & outlet piping
- Fernco® couplers
- Floor / Wall mount bracket
- Services 1 - 10 chairs
- Container Reorder No. **NXT-HG5-002CR**



**NXT Hg5 High Volume
Amalgam Separator**

Item # **NXT-HG5-HV**

Specifications:

- 29¾”h x 17⅝”w x 22½”d
- 3” inlet & outlet piping
- Wall mount bracket
- Services 11 - 20 chairs
- Container Reorder No. **NXT-HG5-002CR**

Hg5 Series of Collection Containers



**NXT Hg5
Collection Container**

with Recycle Kit

Item # **NXT-HG5-002CR**



**Hg5 [Legacy]
Collection Container**

with Recycle Kit

Item # **HG5-002CR**

All Collection Container with Recycle Kits contain:

- Collection Container
- Proper packaging for shipment to our EPA approved recycling facility
- Recycling documentation available online 24/7
- *No contracts!*

Hg5 Series of Compliance Kits



**NXT Hg5
Compliance Kit**

Item # **NXT-HG5-CK**



**Hg5 [Legacy]
Compliance Kit**

Item # **HG5-CK**

All Compliance Kits contain:

- One Collection Container
- One 1.25 gal. Practice Waste Solutions Amalgam Bucket
- One 64 oz. bottle of PowerScrub™ Vacuum Line Cleaner
- Proper packaging for shipment to our EPA approved recycling facility
- Recycling documentation available online 24/7

Solmetex -- Clearly the leader in amalgam separation. The NXT Hg5 amalgam separator is the industry standard and has won more awards than all competitors combined.





The NXT Hg5 Amalgam Separator

Dentistry's Best Selling, Longest Lasting
Amalgam Separator

Trust the Leader

For over 25 years, Solmetex has lead the amalgam separation industry with its award winning NXT Hg5 Amalgam Separator which easily and affordably removes amalgam waste from dental practices.

As a "Total Solution Provider", Solmetex not only takes care of all amalgam waste handling and recycling in the dental practice, but also houses your Certificates of Recycling in our Compliance Center available online 24/7 at solmetex.com.

The NXT Hg5... Good for you and good for the environment.

Solutions

The NXT Hg5 provides the complete solution:

- Available in 3 sizes to suit your practice (mini 1-4 ops; standard 1-10 ops; high volume 11-20 ops)
- Functional with all wet and dry vacuums
- Environmentally friendly packaging
- Simple mail-back recycling program
- Meets all EPA "Dental Rule" requirements
- ISO 11143:2008 certified by NSF International



Easy

- Clear system, easy to inspect
- Two-minute container change out
- No tools necessary

Flexible

- Versatile, wall or floor mount
- Compatible with wet or dry vacuums
- May be installed left or right orientation

Support

- Compliance center
- Online certificates 24/7
- Full Service Customer Support Department

Recycle

- Ship container via UPS
- EPA certified mercury facility
- Recycle Kit includes container and recycle packaging