



# The NXT Hg5<sup>®</sup> Amalgam Separator

## Dentistry's Best Selling, Longest Lasting Amalgam Separator

### Trust the Leader

For over 30 years, Solmetex has helped dental practices safely capture and recycle mercury-laden amalgam with our industry-leading NXT Hg5 Amalgam Separator.

We handle the entire process—from installation to recycling—while keeping your compliance documentation accessible 24/7 through My Solutions Center at [solmetex.com](http://solmetex.com).

### Solutions

#### The NXT Hg5<sup>®</sup> provides the complete solution:

- Available in 3 sizes to suit your practice (mini 1-4 ops; standard 1-10 ops; high volume 11-20 ops)
- Compatible with all wet and dry vacuums
- Environmentally-friendly packaging
- Simple pre-paid mail-back recycling program
- Meets all EPA "Dental Rule" requirements
- ISO 11143:2008 certified by NSF International



***Protecting your practice,  
your patients, and the  
environment from mercury.***

### Easy

- Clear system, easy to inspect
- Two-minute container change out
- No tools necessary

### Flexible

- Versatile, wall or floor mount
- Compatible with wet or dry vacuums
- Install in left or right orientation

### Support

- Online compliance portal
- Access certificates online 24/7
- Full-service customer support

### Recycle

- Ships via UPS with pre-paid shipping label
- EPA-certified mercury recycling facility
- Includes container and safe shipping packaging

# Ordering Information

→ Call Solmetex **877.207.1551**

→ Online: **solmetex.com**

→ If you do not have a dealer, please visit [solmetex.com](http://solmetex.com) and click on "dealer locator" to find a dental dealer near you.

## NXT Hg5® Series of Amalgam Separators

### NXT Hg5® Mini Amalgam Separator

Item #NXT-HG5-MINI



- 25" h x 11" w x 8" d
- ¾" inlet & outlet piping
- Fernco® couplers
- Floor/Wall mount bracket
- Services 1 - 4 chairs
- Container reorder # **NXT-HG5-002CR**

### NXT Hg5® Amalgam Separator

Item #NXT-HG5-001



- 30" h x 11" w x 8" d
- 1½" inlet & outlet piping
- Fernco® couplers
- Floor/Wall mount bracket
- Services 1 - 10 chairs
- Container reorder # **NXT-HG5-002CR**

### NXT Hg5® High Volume Amalgam Separator

Item #NXT-HG5-HV



- 29¾" h x 17" w x 22½" d
- 3" inlet & outlet piping
- Wall mount bracket
- Services 11 - 20 chairs
- Container reorder # **NXT-HG5-002CR**

## Hg5® Series of Collection Containers



### NXT Hg5® Collection Container

with Recycle Kit  
Item #NXT-HG5-002CR



### Hg5® [Legacy] Collection Container

with Recycle Kit  
Item #NXT-G5-002CR

#### All Collection Container with Recycle Kits include:

- Collection Container
- Proper packaging to ship to our EPA-approved recycling facility
- Ready to print, pre-paid UPS shipping
- Recycling documentation available online 24/7

## Hg5® Series of Compliance Kits



### NXT Hg5® Compliance Kit

Item #NXT-HG5-CK



### Hg5® [Legacy] Compliance Kit

Item #HG5-CK

#### All Compliance Kits include:

- One Collection Container
- One 1.25 gal. Practice Waste Solutions Amalgam Bucket
- One 64 oz. bottle of PowerScrub™ Vacuum Line Cleaner
- Proper packaging to ship to our EPA-approved recycling facility
- Recycling documentation available online 24/7



**Solmetex – Clearly the leader in amalgam separation.**

*The NXT Hg5® amalgam separator is the industry standard and has won more awards than all competitors combined.*

