

# Ready for Inspection?

## Here's Your EPA Dental Rule Compliance Checklist

### E.P.A. DENTAL RULE REGULATION



### TOOLS FOR COMPLIANCE

YUP, GOT IT

NOPE, BETTER CALL SOLMETEX!

#### AMALGAM SEPARATOR

Regulation requires the use of an ISO 11143 certified amalgam separator with a separation rate of 95% or greater.

#### NXT Hg5 Amalgam Separator


☐
☐

#### COLLECTION CONTAINER

Regulation requires collection containers be replaced when sediment level reaches the full line or every 12 months, whichever comes first.

**Best Practice:** Don't accidentally get caught with a full container — Always have an extra collection container on hand.

#### Collection Container Recycle Kit


☐
☐

#### AMALGAM WASTE BUCKET

Regulation requires that all 'dry' amalgam waste be collected in a clearly marked amalgam bucket for recycling. Dry waste includes:

- Contact & non-contact amalgam
- Teeth with amalgam
- Chair-side traps
- Vacuum filters
- Spent amalgam capsules

#### Practice Waste Solutions Amalgam Bucket


☐
☐

#### VACUUM LINE CLEANER

Regulation requires the use of only non-chlorine or non-oxidizing disinfectants with a neutral pH  $\geq 6$  and  $\leq 8$ .

#### PowerScrub™ Vacuum Line Cleaner


☐
☐

#### PROOF OF RECYCLING

Regulation requires that dental practices provide amalgam waste recycling records showing name and address of disposal company, plus maintenance and service completed. Records must be kept for at least three years.

#### Certificates of Compliance


☐
☐

## What's the best way to ensure your offices are always in compliance and inspection-ready?

Sign up for the **Solmetex Compliance Program** and your practices will never be without a collection container, amalgam bucket or line cleaner.....keeping them ALWAYS in compliance. Solmetex will automatically send your selected products according to your desired schedule, which you can adjust anytime to suit the needs of your practices.