



STERISIL[®]

Safe, Simple, Sterisil

SHOCK
 TEST
 MAINTAIN


Citrisil[™]

DENTAL WATERLINE MICROBIAL CONTAMINANT CONTROL TABLET




KEY FEATURES

- Cleans and maintains dental unit waterlines to ≤ 10 CFU/ml
- Shock tablets included to kill pre-existing odor causing bacteria
- Eliminates the need to empty bottles and purge lines at the end of each day
- Safe for dental equipment, patients, and clinicians
- Does not adversely affect composite bond strength
- Maintains antimicrobial tubing
- Does not contain oxidizers like iodine or chlorine
- Available for bottles ranging from .70L through 2L bottles




White Tablet



Blue Tablet

48 MAINTENANCE TABLETS



Shock Treatment

2 SHOCK TABLETS INCLUDED WITH EACH BOX

Visit us online at Sterisil.com





STERISIL[®] *Safe, Simple, Sterisil*



PROTECT YOUR PATIENTS AND STAFF AGAINST BACTERIA

Citrisil™ is the only all-in-one shock and daily maintenance dental waterline tablet. Use the orange shock tablet in your lines monthly; allowing it to work overnight to kill pre-existing odor causing bacteria. Use the white or blue maintenance tablet in your bottle daily to continuously clean and maintain dental waterlines. Drop in a new maintenance tablet every time a bottle is refilled. The tablet effervesces and continually stays in the lines to prevent bacteria growth while also being safe for patients and staff.

Part #		White Maintenance Tablets
C50-W	(48) White Maintenance Tablets & (2) Shock Tablets, for .70–1L bottles	
C50-2W	(48) White Maintenance Tablets & (2) Shock Tablets, for bottles ranging from 1L+ to 2L capacity	
Part #		Blue Maintenance Tablets
C50-B	(48) Blue Maintenance Tablets & (2) Shock Tablets, for .70–1L bottles	
C50-2B	(48) White Maintenance Tablets & (2) Shock Tablets, for bottles ranging from 1L+ to 2L capacity	
Part #		Shock Tablets
ECS-20	(20) Shock Tablets	
ECS-4	(4) Shock Tablets	




Proudly made in the USA.

Sterisil, Inc.
835 S. Hwy 105 Ste. D
Palmer Lake, CO USA 80133
+1 719 622 7200

EPA Reg. No. 83315-1
EPA Est. No. 83315-CO-001
Patent No. 6,991,736



1639



EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands



For more information contact your **Sterisil® sales representative or local dealer.**

719 622 7200 or online at Sterisil.com

