



**STERISIL**<sup>®</sup> *Safe, Simple, Sterisil*

○ SHOCK  
 ✓ TEST  
 ○ MAINTAIN

# Dental Unit Waterline Testing

A crucial step in any waterline management program is to perform a water test. Whether you are getting a baseline or validating a waterline treatment regimen, bacteria testing is the only way to know the dental water quality and ensure compliance with the CDC guideline of <500 CFU/ml.

## BACTERIAL ASSESSMENTS FOR ANY SITUATION

### IN-OFFICE TESTS

In-office tests are designed as a screening tool that provides a pass/fail result.



#### RECOMMENDED USE:

- Monthly Testing Interval

PN: MC-4, MC-12



### MAIL-IN TESTS

Mail-in tests are processed by Agencis Labs<sup>™</sup> and include certified results and detailed reports.



#### RECOMMENDED USE:

- Quarterly Testing Interval

PN: D-R2A-1, D-R2A-4, D-R2A-6, D-R2A-10



#### RECOMMENDED USE:

- Quarterly Testing Interval or Failure/Protocol Breach response

PN: D-RC-1, D-RC-4, D-RC-6, D-RC-10



Visit us online at [Sterisil.com](http://Sterisil.com)



# Water Testing Made Easy and Reliable

Pick the test method that matches your chosen testing interval

## MONTHLY



*Easy, In-Office  
Screening Paddle*

- Affordable and proven technology
- No shipping required
- Results in 72hrs
- Easy to use

## QUARTERLY



*The “Gold” Standard  
Mail-in Test*

- EPA validated test method
- Certified 3rd party results
- In-depth reporting with hi-res images of each sample


## QUARTERLY/AS-NEEDED



*Express 24-Hour  
Mail-in Test*

- Results in 24hrs
- A wider spectrum of bacteria captured compared to other test methods
- Certified third-party lab report

The **Organization for Safety, Asepsis, and Prevention** (OSAP) recommends testing monthly until you've attained passing results for 2 consecutive months. Continue to test monthly or reduce to quarterly testing cycles so long as you maintain passing results.



**Now available through your local Sterisil distributor**



For more information contact your **Sterisil<sup>®</sup> sales representative or local dealer.**

719 622 7200 or online at [Sterisil.com](http://Sterisil.com)

