

# SMART Compliance Kits



## Single Operatory Compliance Kit For Bottle-Fed Dental Units



1 Part Number, 3 Products...

*A Complete Waterline  
Compliance Solution*

### EACH SINGLE OPERATORY KIT INCLUDES

- (1) 4-pack Citrisil™ shock tablets
- (1) 4-pack MyCheck™ in-office water test paddles
- (1) Sterisil® Dental Water Microbiological Cartridge, Model Type: Straw
- (1) Quick Start Guide

### KEY FEATURES

- Includes Sterisil's best-selling 365-day straw product w/ built-in shock treatment
- Quick Start Guide condenses the instructions for use (IFUs) for 3 products into one integrated protocol\*
- End-to-end waterline compliance support with one manufacturer
- In-office test paddles provide necessary feedback on the effectiveness of your treatment regimen
- Shock and maintenance products are EPA registered and deliver handpiece water treated to  $\leq 10\text{CFU/ml}$

Visit us online at [Sterisil.com](http://Sterisil.com)

# SAVE \$36.46 INSTANTLY WHEN YOU BUY THE KIT!



**CITRISIL™ SHOCK**  
**SHOCK**

\*\*Price per 4 Tablets:  
**\$10.47\*\***



**MYCHECK™**  
**TEST**

Retail Price:  
**\$88.99**



**STERISIL® STRAW**  
**MAINTAIN**

Retail Price:  
**\$288.99**

RETAIL TOTAL: **\$388.45**

SMART KIT MSRP: **\$351.99**



**SAVE UP TO \$497.81**  
**WHEN YOU BUY THE KIT WITH BUY 3,**  
**GET 1 PRICING OVER PRODUCTS**  
**PURCHASED SEPARATELY AT THE**  
**EQUIVALENT RETAIL PRICE**



Request your **FREE** Waterline Compliance Binder

Contains product literature, instructions for use, safety data sheets, water test logs, and much more!

STERISIL PART #	TYPE	DESCRIPTION
<b>SK-S365</b>	Distilled - SMART Kit	For bottle-fed delivery units using distilled water or tap water containing 0-100ppm TDS*. Includes Citrisil™ shock, 365-day straw, and MyCheck™ in-office water test paddles.
<b>SK-S365M</b>	Municipal - SMART Kit	For bottle-fed delivery units using municipal tap water containing 101-250ppm TDS*. Includes Citrisil™ shock, 365-day straw, and MyCheck™ in-office water test paddles.



Proudly made  
in the USA.

Sterisil, Inc.  
835 S. Hwy 105 Ste. D  
Palmer Lake, CO USA 80133  
+1 719 622 7200  
EPA Reg. No. 83315-1  
EPA Reg. No. 83315-2  
Patent No. 6,991,736

\*If you intend to use municipal tap water in your dental unit, we recommend assessing the TDS content prior to purchasing your kit. Clinics with >250ppm TDS should use distilled water and PN: SK-S365.



For more information contact your **Sterisil® sales representative or local dealer.**

719 622 7200 or online at [Sterisil.com](http://Sterisil.com)