



SMART

Compliance Kits

Supplemental Shock & Test Kits

EASILY ADD WATER TESTING AND SHOCK TO EXISTING MAINTENANCE REGIMENS



- Each supplemental test kit includes:
 - (1) 4-pack Citrisil™ shock tablets
 - (1) 4-pack MyCheck™ water test paddles
 - (1) 4-vial mail-in water test kit
- Shock and testing guide makes integrating these protocols easy and seamless with any existing regimen
- End-to-end waterline compliance support with one manufacturer
- Non-oxidizing shock treatment
- Citrisil™ shock tablets are EPA registered and designed to kill odor causing bacteria that contribute to water test failures
- Choose from 2 different kit options based on the mail-in test method used:
 - Sterisil **SMART** 24HR Water Testing Kit: SK-RC-4
 - Sterisil **SMART** R2A Water Testing Kit: SK-R2A-4



Visit us online at Sterisil.com



Complete the Cycle by Adding Shock and Testing to Your Existing Treatment Regimen

SHOCK TREATMENTS ARE IMPORTANT. HERE'S WHY:

- Extreme Variability in Source Water Quality**
 Erodes the effectiveness of maintenance treatments and allows bacteria to recolonize waterlines.
- Human Error**
 Common in protocol driven routines. Monthly or quarterly shock cycles reduce the impact of undetected protocol breaches that would otherwise lead to a test failure.
- Water Test Failures**
 Your response will always be the same, shock immediately and retest ASAP.


WATER TESTING IS IMPORTANT. HERE'S WHY:

- Not Testing is Not Knowing**
 If you don't assess the bacterial content in your dental unit, how do you know you're compliant?
- Liability Claims**
 Patient injury caused by questionable infection control practices can open the door for lawsuits*. Consistent and properly documented water testing provides the basis for a legal defense against claims of negligence should the quality of your dental unit water be called into question.
- Legal Requirements**
 In some states, clinicians are now required to prove the bacterial content in dental unit waterlines meets the <500 CFU/ml drinking water standard.

STERISIL PART #	TYPE	DESCRIPTION
SK-RC-4	Sterisil SMART 24HR Water Testing Kit	SMART Kit intended to supplement water testing and shock for any bottle-fed treatment regimen. Includes (1) 4-vial RapidCheck™ test kit and (1) 4-pack of Citrisil™ shock tablets. Maintenance treatment not included.
SK-R2A-4	Sterisil SMART R2A Water Testing Kit	SMART Kit intended to supplement water testing and shock for any bottle-fed treatment regimen. Includes (1) 4-vial StandardCheck™ test kit and (1) 4-pack of Citrisil™ shock tablets. Maintenance treatment not included.



Proudly made in the USA.

 Sterisil, Inc.
 835 S. Hwy 105 Ste. D
 Palmer Lake, CO USA 80133
 +1 719 622 7200

EPA Reg. No. 83315-1
 Patent No. 6,991,736

*<https://www.rdhmag.com/pathology/oral-pathology/article/16406303/top-reasons-hygienists-are-sued>



For more information contact your **Sterisil® sales representative or local dealer.**

719 622 7200 or online at Sterisil.com

