



BEFORE.



AFTER.

Isolation. **Simplified.**



DryShield®



DISCOVER THE DRYSHIELD DIFFERENCE

Practice dentistry like never before.

High suction. Saliva ejectors. Bite blocks. Dry angles. Cotton rolls. Gauze. Replace them all with one simple, powerful isolation solution. The patented DryShield system was designed with the dentist in mind, offering total control of your patient's mouth by keeping the working area dry, clean and visible. Remarkably, DryShield dentists are seeing chair time reduced by an average of 25-40% per procedure. With more free time and less variability, doctors and hygienists are doing more procedures in each appointment, scheduling more patients in a day, and consistently running their schedules on time. At the same time, patients are enjoying a more comfortable experience as well as shorter appointments.

ONE SIMPLE DEVICE THAT DOES IT ALL

DryShield performs all of the tasks of a high-suction evacuator, bite block, tongue guard and oral pathway protector in one easy-to-use device. It's an intelligent breakthrough in isolation, bringing unmatched levels of ease and comfort to doctors and patients alike. Opt for the cost savings of autoclavable mouthpieces, or the convenience of single-use mouthpieces.

Autoclavable Mouthpieces



- 1 Cheek Retractor**
Keeps patient's cheek retracted to improve visibility and protect it from trauma
- 2 Airway Protection**
Protects airway from particles and debris
- 3 Tongue and Oral Pathway Shield**
Keeps patient's tongue shielded and oral pathway protected from trauma
- 4 Bite Block**
Autoclavable: Interchangeable bite blocks available in multiple sizes to facilitate access and maximize patient comfort for all mouth sizes
Single-Use: Built-in bite blocks maximize patient comfort
- 5 Suction**
Removes water, saliva and debris from patient's mouth to keep your work area dry
- 6 Easy-to-Insert Design**
Mouthpiece can be inserted in a few seconds, and removed even faster

Single-Use Mouthpieces



SO MUCH TO LOVE, SO EASY TO USE

What sets DryShield apart?

- DryShield has engineered two options — the Autoclavable and Single-Use mouthpieces — to meet your needs, for maximum affordability or convenience.
- Installation is simple, attaching directly to the HVE — no special equipment required.
- All DryShield mouthpieces are made from soft and flexible material that contours to the patient's mouth.
- Each mouthpiece size features a matching bite block, giving the patient a firm place to rest for a perfect fit in every mouth.
- Hands-free, continuous suction allows patients to be comfortable and assistants to multi-task.
- The DryShield system is portable between operatories and offices.

SEAMLESS, COST-EFFECTIVE INTEGRATION

Fully portable. Installs in seconds.



Find out more at DryShield.com or 1 888 DRY 0300



11230 Grace Avenue, Suite B | Fountain Valley, California 92708 | 1 888 379 0300 | DryShield.com

DryShield®

© 2019 DryShield, LLC. All rights reserved. DryShield and the DryShield logo are trademarks of DryShield, LLC. Patented and patent pending in the US and worldwide. DryShield.com/patents