




STERISIL®

Safe, Simple, Sterisil

✓ **SHOCK**
○ **TEST**
○ **MAINTAIN**

Citrisil™ Shock Tablets

ELIMINATING WATERLINE CONTAMINATION,
ONE TABLET AT A TIME



PNs: ECS-20

KEY FEATURES

- Cleans dental unit waterlines
- Kills odor causing bacteria
- EPA claim: ≤ 10 CFU/ml
- Cleans the inner surface of dental bottles
- Orange colorant makes the shock visually apparent in the bottle
- Non-oxidizing formula
- Chemically compatible with the Sterisil® Antimicrobial Bottle With Built-In Antimicrobial Properties To Protect The Bottle

Visit us online at Sterisil.com



PROTECT YOUR COMPLIANCE STATUS FROM COMMON PITFALLS:

PROTOCOL BREACHES

Any number of off label practices can put waterline compliance in jeopardy.

WATER TEST FAILURES

Regardless of why you failed, the response will always be the same; shock immediately.


INCONSISTENT SOURCE WATER QUALITY

If you use municipal water, you have more impurities and microbes entering the system that feed biofilm development.

PLUMBING OR EQUIPMENT RELATED EVENTS

Regular occurrences like boil water advisories or routine chair maintenance may necessitate a shock.

TAKE CHARGE OF YOUR BOTTLE HYGIENE



Use Citrisil™ Shock tablets to clean the bottle's inner surface every time you shock your dental unit waterlines.



STERISIL PART #	DESCRIPTION
ECS-20	(20) Citrisil™ Shock tablets




Proudly made
in the USA.

Sterisil, Inc.
835 S. Hwy 105 Ste. D
Palmer Lake, CO USA 80133
+1 719 622 7200

EPA Reg. No. 83315-1
EPA Est. No. 83315-CO-001
Patent No. 6,991,736



1639



EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands



For more information contact your **Sterisil® sales representative or local dealer.**
719 622 7200 or online at Sterisil.com

