



STERISIL[®] *Safe, Simple, Sterisil*

**AUTOCLAVE
DENTAL**

Sterisil[®] System G4

AUTOCLAVE AND DENTAL WATER,
MADE EASY AND RELIABLE



KEY FEATURES

- Effluent dental water treated to $\leq 10\text{CFU/ml}$
- Equipped with a Class B UV light
- Eliminates the need for purchasing and storing distilled water
- Output autoclave water is purified to $<10\text{ppm TDS}$ and is suitable for any washer or autoclave
- Live water quality monitoring
- Audible and visual alarms
- Self-closing valve heads make routine maintenance easy
- User serviceability eliminates costly technician visits for routine maintenance

The Sterisil[®] System G4 is designed to deliver purified autoclave water as well as treated dental water, made right from your municipal supply. A simple plug and play solution to all of dentistry's water related supply problems.

Visit us online at Sterisil.com

