

DryShield® All-in-One Isolation System

ISOLATION • PROTECTION • WORKFLOW • SIMPLIFIED



DryShield®



BENEFITS

- Cut 25-40% in procedure time
- More production AND higher quality of work
- Two-quadrant dentistry increases visibility and reduces chair time
- Protect patients and staff from infectious aerosols
- Clinicians can work more independently
- Hands-free operation allows patients to be comfortable and assistants to multi-task

COMPONENTS

- Cheek retractor improves visibility and protects soft tissue from trauma
- Airway protection prevents aspiration of foreign objects
- Tongue and oral pathway shield keeps tongue shielded and oral pathway protected from trauma
- Bite block maximizes patient comfort (Autoclavable bite blocks are interchangeable to facilitate further size customization)
- Constant suction keeps the work area dry and removes aerosol right at its source
- Easy-to-insert design can be installed and removed from patient's mouth in seconds
- Made from soft and flexible material that contours to the patient's mouth
- System is portable between operatories and offices

Discover the DryShield Difference

DryShield performs all of the tasks of a high-suction evacuator, bite block, tongue guard and oral pathway protector in one easy-to-use device. The patented DryShield system was designed with the dentist in mind, offering total control of your patient's mouth by keeping the working area dry, clean and visible.

It's an intelligent breakthrough in isolation, bringing unmatched levels of ease and comfort to doctors and patients alike. Opt for the cost savings of autoclavable mouthpieces, or the convenience of single-use mouthpieces.

Autoclavable and Single-Use Mouthpiece Options:



- Autoclavable Mouthpieces
- 4 mouthpieces per box
- Autoclavable up to 50 times each
- Removable bite blocks for custom sizing



- Single-Use Mouthpieces
- 20 mouthpieces per box
- Disposable mouthpiece option
- Built-in bite blocks maximize patient comfort

Visit us online at DryShield.com



SAFE WATER



SAFE PATIENT



SAFE ENVIRONMENT

Dentistry's Leading Provider of Water Safety Solutions

PART #	PRODUCT NAME	DESCRIPTION
DS-SYS-001	DryShield Isolation System	All-in-one isolation system, Includes: DS1 Isolation System, Portable Operatory Kit, Y Connector Kit, Cleaning Brush, and additional O-Ring
DS-SYS-002	DS Lite Adapter	DryShield DS1 adapter
DS-MPL-005	DryShield Autoclavable Mouthpiece - Large	Large autoclavable silicone mouthpieces; ideal for adults with a large vertical opening
DS-MPM-004	DryShield Autoclavable Mouthpiece - Medium	Medium autoclavable silicone mouthpieces; ideal for adults or teens with existing permanent 2nd molars
DS-MPS-003	DryShield Autoclavable Mouthpiece - Small	Small autoclavable silicone mouthpieces; ideal for children 9 years of age to young teens/adult with full-erupted permanent 1st molars
DS-MPP-002	DryShield Autoclavable Mouthpiece - Pedo	Pedo size autoclavable silicone mouthpieces; ideal for young children 8 years of age or younger, with primary or mixed dentition
DS-MPX-001	DryShield Autoclavable Mouthpiece - X-Pedo	X-Pedo size autoclavable silicone mouthpieces; ideal for very young children 3-4 years or younger
DS-SUL-600	DryShield Single-Use Mouthpiece - Large	Large single-use mouthpieces; ideal for adults with a large vertical opening
DS-SUM-500	DryShield Single-Use Mouthpiece - Medium	Medium single-use mouthpieces; ideal for adults or teens with existing permanent 2nd molars
DS-SUS-400	DryShield Single-Use Mouthpiece - Small	Small single-use mouthpieces; ideal for children 9 years of age to young teens/adult with full-erupted permanent 1st molars
DS-SUP-300	DryShield Single-Use Mouthpiece - Pedo	Pedo size single-use mouthpieces; ideal for young children 8 years of age or younger, with primary or mixed dentition
DS-SK-001	DryShield Starter Kit	Everything You Need to Get Started Includes: DS1 Isolation System, Portable Operatory Kit, Y Connector Kit, Cleaning Brush, and 4 boxes of autoclavable mouthpieces (Lg, Md, Sm, Pd)"

Visit us online at DryShield.com