



SHOCK     TEST     MAINTAIN

# Sterisil® R2A Waterline Test Kit

## Mail-In Water Testing The Gold Standard R2A Test



### BENEFITS

- EPA certified third-party lab results
- Post result consultation
- Results in 5–7 days

### FEATURES

- R2A Heterotrophic Plate Count (HPC) Water Test
- EPA method spread plating
- Automatic digital counting and recording
- Results posted in My Solutions Center

Solmetex  
System of  
Solutions

Source  
Water  
Quality



Dental Unit  
Waterline  
Management



Isolation  
and Aerosol  
Control



Evacuation



Amalgam  
Capture and  
Recycling



Waste  
Management



Every  
Drop  
Matters™





Sterisil R2A bacterial tests utilizes an EPA validated R2A plating method to grow/count HPC bacterial colonies in samples collected from dental unit waterlines. Tests are conducted at EPA Certified lab, Dentisafe (37917). Kit includes vials, ice packs, return packaging, chain of custody and instructions, next day air return shipping label. Results in days.



**TESTING IS SIMPLE! SCAN HERE FOR INSTRUCTIONS**  
[Bit.ly/fastcheck15-ifu](https://bit.ly/fastcheck15-ifu)



**It is important to document your test results.**  
Sterisil recommends that you log them at My Solutions Center.  
[Bit.ly/mysolutionscenter](https://bit.ly/mysolutionscenter)



**Scan to schedule a Free Water Consultation with one of our experienced Water Specialists.**  
Failed the test? Chat with us. Scan the code.  
[Bit.ly/water-consult](https://bit.ly/water-consult)

STERISIL PART #	SHELF LIFE	DESCRIPTION
SDENT-R2A-1	Unlimited	VIAL COUNT 1
SDENT-R2A-4	Unlimited	VIAL COUNT 4
SDENT-R2A-6	Unlimited	VIAL COUNT 6
SDENT-R2A-8	Unlimited	VIAL COUNT 8
SDENT-R2A-12	Unlimited	VIAL COUNT 12
SDENT-R2A-16	Unlimited	VIAL COUNT 16

**Order/Info:**

**(800) 216-5505**  
**Solmetex.com**



**Solmetex**  
**Sterisil**  
**DryShield**